

ESC derivation – Paul updated 1-5-17

- *I combined different protocols (see supplemental protocols). The major changes are:*
 - a. *Using 2i/LIF medium (M15 medium, 1 μ M PD0325901 [Sigma #PZ0162], 3 μ M CHIR99021 [Millipore #361559], 100U/ml LIF [1/10 LIF used for regular ES culture]) and MEF feeders during the first few (P1-P3) passages*
 - b. *Using accutase (Thermo Cat. # A1110501) to dissociate blastocysts*
1. Isolate blastocysts from uteri following “Manipulating the Mouse Embryo” textbook (check the supplemental protocols for more detail)
 - a. Before opening the mouse, put 2i/L on top of MEF feeders on 12-well plate (I usually prepare 8 wells). Pre-warm in the CO₂ incubator for at least 30 min.
 - b. Put 5ml Knockout DMEM (Thermo #10829018) + 20 mM HEPES on each of the three 6 cm plates (Thermo #15630106). Pre-warm in the CO₂ incubator for at least 30 min before opening the mouse.
 - c. Isolate uteri from the mouse. Cut off oviducts and ovaries from both ends. Then cut the “Y”-shaped uterus from the middle to break it into two uterine horns.
 - d. Flush the uterine horns using syringe from both ends with the pre-warmed Knockout DMEM+HEPES (use 4-8 ml in total).
 - e. Transfer individual blastocyst to each well of the 12-well plate (with MEF feeders and 2i/L). If there are >8 embryos, put 2 blastocysts in one well.
 2. Culture blastocysts in 2i/L without ANY disturbance for 4 days. After the 4th day, check the outgrowth under the microscope. The blastocysts with outgrowth can be trypsinized and transferred at 5-6th days following **Step 3**.
 - a. Blastocysts hatched faster need to be transferred earlier.
 - b. To find the blastocyst, look through the bottom of the plate towards the light. The blastocysts appear as white dots and can be circled with the pen.
 - c. 80-90% blastocysts should give rise to outgrowth.
 3. Handpick blastocysts with outgrowth (P1)

- a. Pre-warm ESC medium (M15+LIF) on top of feeders growing on the 24-well plate. The number of wells = the number of blastocysts to be passaged.
 - b. Don't change the medium (2i/L) before handpicking. The blastocysts are loosely attached to the feeder. Also avoid shaking while moving the plate.
 - c. Put 30 μ l accutase (don't warm up, directly take from refrigerator and bring it to the room temperature. Warming up accutase decreases its activity) into eppendorfs (= the number of blastocysts going to be picked)
 - d. Handpick the colonies with 10 μ l tip. Directly transfer the blastocyst to the accutase in the eppendorf and pipet several times to ensure the transfer. It's better to pick all the colonies asap (within 5-10 min).
 - e. Incubate all the eppendorfs in the CO₂ incubator for 15-20 min.
 - f. Use 100 μ l tip to put 100 μ l M15 medium to the eppendorf, pipet up and down 5-6 times until blastocysts are dissociated. Don't pipette more than 10 times. Cells still aggregate together can be dissociated at P2.
 - g. Centrifuge at 800g for 5 min
 - h. Carefully remove the top 100 μ l medium from the eppendorf. Use the pre-warmed ESC medium on 24-well plate (**Step 3a**) to resuspend the pellet and put back into the same well.
 - i. Rock the plate to evenly distribute cells. Grow cells in the CO₂ incubator.
4. Change the ESC medium daily until the colony are visible with appropriate size
 - a. Blastocysts that not dissociated well (**Step 3f**) can be trypsinized with 0.25% TE (red TE) and passaged to the new 24-well plate the 1-2 days after P1
 5. Trypsinize and transfer (using 0.25% TE) the blastocysts on 24-well plate until it reaches 50-60% confluency. Save ESC stocks at early passages. Then ESC lines can be expanded for experiments.