# <u>General IHC Protocol for Formalin-fixed Paraffin Embedded Staining</u> (Steps 16-18: follow only if using biotinylated secondary Ab)

### **Materials**

- 1. Histoclear
- 2. Ethanol (100%, 95%, 70%)
- 3. Albumin from bovine serum (BSA): Sigma A7906-500g
- 4. Primary Antibody
- 5. Anti-primary (Secondary) Antibody
- 6. Permount mounting medium
- 7. Triton X-100: Sigma T8787-100ml

- 8. Modified Meyer's Hematoxylin
- 9. Phosphate Buffered Saline 1x
- 10.  $H_2O_2$  30% solution
- 11. Coplin jars (plastic for Ag retrieval)
- 12. Pressure cooker
- 13. Dako pen
- 14. DAB substrate reagents
- 15. Slide Coverslips

# **Solutions**

- 1. 5% BSA/0.3% Triton X-100 (should be made fresh every 1-2 months)
  - a. 25 g of BSA into 500 ml of 1x PBS
  - b. Add 1.5 ml of 100% Triton X-100
  - c. Filter and store at 4°C
- 2. 1% BSA/0.3% Triton (should be made fresh every 1-2 months)
  - a. 5 g of BSA into 500 ml of 1x PBS
  - b. Add 1.5 ml of 100% Triton X-100
  - c. Filter and store at 4°C
- 3. Sodium Citrate pH 6.0 (Genentech Media Prep), or DAKO Antigen Retrieval Buffer (see step 6a)

# **Procedure**

- 1. Fix tissues in formalin immediately after collection 1hr/mm of tissue (perfused tissue is best)
- 2. Wash in PBS 1x for 4x15min and store in 70% EtOH
- 3. Embed in paraffin and section

# <u>Day 1 (~6 hours)</u>

5. Deparaffinize slides 60min:

Done using 250 ml of solution in staining buckets with slide rack

Note – never let your slides dry out!!!!!! Throughout this whole procedure



- 6. Antigen retrieval (Heat Induced Epitope Retrieval) 30min:
  - a. In Plastic Coplin jar(s), place slides in <u>Dako Antigen Retrieval Buffer</u> (To prep: in fridge, make excess of 1:10 dilution with distilled water, pH to 6.0-6.2, optimal is pH 6.1, using HCl if necessary.) (Alternatively, may use 10 mM <u>sodium citrate buffer</u>); antigen retrieve for 10 min at low pressure in the pressure cooker (filled between the 6 and 8 mark with MQ water).
    - i. When pressure cooker is done (*takes roughly 20 minutes*) quick release the steam, allowing pressure to normalize
    - ii. Take off lid, remove Coplin jar(s) careful! They'll be hot, use oven mitt and set them on bench, remove lids and allow them to cool to ~tepid RT. Wait for water to cool before discarding from the cooker.
    - iii. Pull Coplin jar out and let cool to warmish RT 20min maybe?
- 7. Wash in PBS 1x5 min in Coplin jar, on bench, RT
- 8. Incubate in 3% H2O2 10 min RT (25mls 30% H2O2 up to 250mls MQH20) in Coplin jar, on bench
- 9. Wash with PBS 5 min in Coplin jar, on bench
- 10. Block for 1-3 hours at room temperature 5% BSA/0.3% Triton X-100 (If having background issues, 3-6hrs works well). Block in the Coplin jars either on bench top or (if background issues) on the Bellydance shaker
- 11. Primary antibodies overnight (>10 hrs) at 4°C in 1% BSA/0.3% Triton X-100. Use 200ul of primary antibody solution per slide (if smaller tissues, use ~50ul). Best if in a level, moist container (blue slide containers, holds 20 slides each, or slide boxes with 1" MQ water in bottom).

| Primary Antibody | Concentration |
|------------------|---------------|
| Casp3            | 1:400         |
| рН3              | 1:300         |
| Pax 5            | 1:200         |

### \*\*Do NOT add Ab to your unstained control here!

Instead, incubate in 1% BSA/0.3% Triton X-100 o/n

### Day 2(~5.5 hrs-6.5 hrs)

- 12. Save primary Abs label with concentration and date (store at 4C for 2wks-1mo)
- 13. Rinse in PBS 4 x 15 min in Coplin jars on Belly-dance shaker at room temperature (60min)
- 14. Secondary antibodies (1:500) for 2 hrs at room temperature in 1% BSA/0.3% Triton X-100. [250ul (or~50ul) per slide] Use blue slide containers or slide boxes.

| Primary Antibody | Secondary Antibody<br>(1:500) |
|------------------|-------------------------------|
| Casp3            | Biotinylated anti-rabbit      |
| pH3              | Anti-rabbit                   |
| Pax 5            | Biotinylated anti-goat        |

15. Wash in PBS 4 x 15min (60 min) in Coplin jar on Belly-dance shaker.

#### **STEPS 16-18 ONLY APPLY TO RXNS WITH BIOTIN**

- 16. Secondary antibodies Biotin conjugated 2° AB
- 17. Incubate for 30min in ABC reagent (use blue slide containers/ slide boxes)
- 18. Wash for 2x5min in PBS

#### 19. Perform DAB reaction

- Make DAB reagent (400ul per slide)
- Make sure you have a bucket filled with MilliQ  $H_20$  ready
  - \*\*Cover 15ml falcon tube in foil reagent needs to be kept in dark and used within 30min\*\*
  - 1 drop Vector A per 1ml MQ H20 mix
  - 2 drops Vector B per 1 ml
  - 1 drop Vector C per 1 ml
  - Mix by inversion, wait 2 mins before adding to samples.

Incubate slides for 10 min sharp – in the dark (foil) – used inverted pipet tip boxes to keep in the dark.

Usually do slides 2 by 2 so that it is more precisely 10 min

- 20. Place immediately in the H2O after 5min of DAB leave there until all slides have finished the DAB incubation
- 21. Change H<sub>2</sub>O to fresh H<sub>2</sub>O wash 5min if slides look particularly brown do this again
- 22. Counterstain with modified Meyer's Hematoxylin 10 min

Put in warm tap water bucket and run warm tap water over it for 8 min

#### 23. Dehydrate:

70%EtOH x 3min 95%EtOH x 3min 95%EtOH x 3min 100%EtOH x 3min 100% EtOH x 3min Histoclear x 10min Histoclear x 10min

Coverslip with Permount Use about 100ul Permount per slide Wipe off excess Let dry at least 1 hour before looking at on scope