K-ras genotyping (LSL, G12D)

Note: Carl has done seperate pcr using U1-wt; and U1-Mut as sepearte PCR. 3 primer PCR does not work well with Carl and Marcelo. Although it has worked well with Jacks lab and Lowe lab. One should always try the three primer system before using the two primer system.

Notes:

- This reaction looks for the presence of the stop cassette NOT for the mutation. Samples must be sequenced to look for mutation.
- Run a large gel at 100V for at least 40-50 min to allow the two bands to separate well!
- Run 2 % gel

Primers:

K-ras-U1	5' – CGCAGACTGTAGAGCAGCG – 3'
K-ras-wt	5' – GTCGACAAGCTCATGCGGG – 3'
K-ras-mut	5' – CCATGGCTTGAGTAAGTCTGC – 3'

Size of bands:

 Wt
 500 bp

 Mut
 550 bp

PCR conditions:

94 °C for 3min 94 °C for 30 sec 60 °C for 1:30min 72 °C for 1min 30 cycles 72 °C for 5min 10C for ever

50µl reaction

2µl DNA 1µl K-ras-U1 1µl K-ras-wt 1µl K-ras-mut 1µl dNTPs (10mM) 1µl Taq 5µl 10x buffer add H2O

20ul reaction (Bin Xue using GoTaq to give more reliable results)

1 μl DNA 0.5 μl K-ras-U1 0.5 μl K-ras-mut 10 μl GoTaq 8 µl H2O