## Osteoclast Differentiation (in vitro)

## Day 1:

Sacrifice & dissect:

- 1. Sacrifice mouse
- 2. Dissect out both femurs and tibias

3. Dip bone briefly in 70% EtOH then place in sterile PBS in 6 well plate

Isolate cells:

- 1. Prefill 10 mL syringe with 25G needle with cold sterile PBS.
- 2. In sterile 6-well plate, cut off both ends of bone to be flushed. For tibia cut off right where you can see the red marrow.
- 3. Hold bone with forceps and insert 25G needle into bone end, crunch up and down while flushing PBS through bone into an empty well of 6-well plate. Flip bone over and repeat on other end. Repeat procedure until bone is white. If cut off ends look like they have marrow, flush them also.
- 4. Pipet PBS + marrow repeatedly and transfer to a 50 mL conical.
- 5. Spin cells at 1000rpm for 5 minute
- 6. Add 5 mL RBC lysis buffer, incubate RT for ~5 minutes, then spin at 1000rpm for 5 minutes.
- 9. Resuspend cells in α-MEM media complete + MCSF (10ng/ml) or 5% CMG; Plate cells at 2 legs/ 6well plates with 6-8 mL media per well

Day 2: Harvest Non-Adherent cells (bone marrow macrophages):

- 1. Stromal cells and some macrophages will adhere to the plastic within 4 hours. You can collect non-adherent cells for differentiation into macrophages or OCL anytime between 4hrs to day 3. The longer you wait, the lower your yield of nonadherent cells in my experience. I usually plate bone marrow in the evening and collect non-adherent cells after overnight culture first thing in the morning (12-15 hrs).
- 2. Collect non-adherent cells and count
- 3. Spin cells and resuspend in  $\alpha$ -MEM media + MCSF (10 ng/ml or 5% CMG).
  - a. Plate **125K cells**/well of 48 well plates
    - b. Volume of media should be 200ul/96well; 0.5ml/48 well; 2ml/24 well; 4-5ml/6well
    - c. pH of media should be 7.2-7.4 for good osteoclast development
    - d. Plate in multiples so that 2-3 wells are for TRAP staining alone, 1-3 wells for dentine, 1-2 wells for BD biocoat, 3 wells for RNA.

Day 5: Change to fresh media

1. Change media to fresh  $\alpha$ -MEM + MCSF/CMG.

Day 8: Change to fresh media

- 1. Change media to fresh  $\alpha$ -MEM + MCSF/CMG.
- 2. Note: you may change to add RANKL depending on the density of the expanding macrophages; if dense enough you can skip to day 11.

Day 11: Change to fresh media

- 1. Change media to fresh α-MEM + MCSF/CMG + 50-100 ng/mL RANKL.
- 2. Begin to monitor daily for formation of osteoclasts; once large osteoclasts have formed, they will apoptose in 24-48 hours, so fix once large osteoclasts are observed

Day 14: Change to fresh media

- 1. Change media to fresh  $\alpha$ -MEM + MCSF/CMG + 50-100 ng/mL RANKL.
- 2. Continue monitoring for osteoclasts.

## **Fixation for TRAP staining:**

- 1. Remove media and rinse 1x in PBS
- 2. Fix in 3.7% formaldehyde in PBS for 10 minutes at RT.
- 3. Rinse 2x in PBS. Flick to empty plates, **do not** aspirate as this will suck up osteoclasts
- 4. Plates can be stored in PBS/1% sodium azide before staining (if so, remember to rinse with PBS 2x before starting staining as sodium azide inhibits TRAP staining)

## TRAP staining in plates:

- 1. Make up TRAP staining solution as per instructions
- 2. Prewarm plates with PBS at 37°C
- 3. Add substrate to fixed cells and incubate 15 min 1 hr (usually an hour).
- 4. Rinse with water or PBS several times; can be stored in PBS/1% sodium azide at 4°C.
- 5.

However, If good osteoclasts are visible, then consider stopping assays.

- a) TRAP staining see below
- b) resorption assays either change the media to acidified alpha MEM complete +RANKL and MCSF or allow assays to continue without changing media which will acidify naturally. Acidify alpha MEM complete by: aliquot out 50cc to conical tube. Add 1 drop 12N HCL. Check that pH is 6.8-7.0 with pH paper. Sterile filter.