

TRAP Staining for Formalin-fixed Paraffin Embedded Sections

Procedure

1. Pre-warm MQ H₂O (~50ml per coplin jar) at 37-42°C (for step 3b) ~30min.

2. Deparaffinize slides ~60min:

Best if you can heat, but doesn't really matter

Note – never let your slides dry out!!!!!! Throughout this whole procedure

Deparaffinize

1.) Histoclear x 10min

2.) Histoclear x 10min

Rehydrate

3.) 100% EtOH (or Flex 100) x 3 min

4.) 100% EtOH (or Flex 100) x 3 min

5.) 95% EtOH (or Flex 95) x 3 min

6.) 95% EtOH (or Flex 95) x 3 min

7.) 70% EtOH x 5min

8.) MQ H₂O x 5 min

Note all of the above steps do not have to be exact... you can leave for longer

3. Prepare the staining solution (according to Beaker B in steps 3 and 4 of the Acid Phosphatase, Leukocyte Kit).

a. Make 1 mL Diazotized Fast Garnet GBC Solution (1:1 Fast Garnet GBC Base Solution: Sodium Nitrite Solution). Mix by gentle inversion, let stand for 2 min

b. Mix the following solutions while stirring (following volumes = 1 coplin jar):

MQ H₂O pre-warmed to 37°C – **45mL**

Diazotized Fast Garnet GBC Solution (6a) – **1mL**

Naphthol AS-BI Phosphate Solution – **0.5mL**

Acetate Solution – **2mL**

Tartrate Solution – **1mL**

***scale volumes as necessary for container**

c. Pour solution into coplin jar and ensure that temperature is **37-42°C** (Important!)

4. Incubate slides at 37-42°C in staining solution for 1 hour.

5. Rinse slides thoroughly in MQ H₂O (2x 5-10 min).

6. Counterstain with Hematoxylin Solution, Gill No.3 for 2 min.

Put in warm tap water bucket and run warm tap water over it for 8 min

Coverslip with Fluoromount G (Aqueous mounting medium)

Use about 100ul permount per slide

Wipe off excess

Seal with nail polish

Let dry at least 1 hour before looking at on scope (usually O/N)